

# Q-interactive® System Requirements



Q-interactive requires two Apple iPad® devices that use iPadOS 13 or higher and have at least 16GB of storage.

## Tablet Requirements

Family	Model	Practitioner Device		Client Device	
		Supported	Not Supported	Supported	Not Supported
iPad	iPad 1 <sup>st</sup> to 4 <sup>th</sup> gen		X		X
	iPad 5 <sup>th</sup> gen (2017)	✓		✓	
	iPad 6 <sup>th</sup> gen (2018)	✓		✓	
	iPad 7 <sup>th</sup> gen (2019) –10.2” screen	✓		✓	
	iPad 8 <sup>th</sup> gen (2020) –10.2” screen	✓		✓	
	iPad 9 <sup>th</sup> gen (2021) – 10.2” screen	✓		✓	
iPad Air	iPad Air 1 <sup>st</sup> gen		X		X
	iPad Air 2 <sup>nd</sup> gen	✓		✓	
	iPad Air 3 <sup>rd</sup> gen (2019) –10.5” screen	✓		✓	
	iPad Air 4 <sup>th</sup> gen (2020) – 10.9” screen	✓		✓	
iPad mini	iPad mini 1 <sup>st</sup> to 3 <sup>rd</sup> gen		X		X
	iPad mini 4 <sup>th</sup> gen	✓			X*
	iPad mini 5 <sup>th</sup> gen (2019)	✓			X*
	iPad mini 6 <sup>th</sup> gen (2021)	✓			X*
iPad Pro	iPad Pro 9.7”	✓		✓	
	iPad Pro 10.5”	✓		✓	
	iPad Pro 11” (1st, 2nd & 3rd gen)	✓		✓	
	iPad Pro 12.9” (1st to 5th gen)	✓		✓	

\* Not supported because the psychometric impacts of displaying stimuli on a larger or smaller screen are unknown.